

7th Grade Science DNA Extraction Lab

Name _____ Date _____

DNA Extraction: Strawberry

Background: The long, thick fibers of DNA store the information for the functioning of the chemistry of life. DNA is present in every cell of plants and animals. The DNA found in strawberry cells can be extracted using common, everyday materials. We will use an extraction buffer containing salt, to break up protein chains that bind around the nucleic acids, and dish soap to dissolve the lipid (fat) part of the strawberry cell wall and nuclear membrane. This extraction buffer will help provide us access to the DNA inside the cells.

Pre-lab questions:


1. What do you think the DNA will look like?
2. Where is DNA found?

Materials:

heavy duty ziploc bag
1 strawberry
10 mL DNA extraction buffer (soapy, salty water)
cheesecloth
funnel
50mL vial / test tube
glass rod, inoculating loop, or popsicle stick
10 mL ethanol

Procedure:

1. Place one strawberry in a Ziploc bag.
2. Smash/grind up the strawberry using your fist and fingers for 2 minutes. *Careful not to break the bag!!*
3. Add the provided 10mL of extraction buffer (salt and soap solution) to the bag.
4. Knead/mush the strawberry in the bag again for 1 minute very slowly so no suds are seen.
5. Assemble your filtration apparatus as shown to the right. Use a rubber band to secure the cheesecloth.
6. Pour the strawberry slurry into the filtration apparatus and let it drip directly into your test tube.
7. Tilt the test tube and slowly pour 10 ml cold ethanol into the tube. OBSERVE _
8. Dip the loop or glass rod into the tube where the strawberry extract and ethanol layers come into contact with each other. OBSERVE _



Conclusions and Analysis

1. It is important that you understand the steps in the extraction procedure and why each step was necessary. Each step in the procedure aided in isolating the DNA from other cellular

materials. Match the procedure with its function:

PROCEDURE FUNCTION

- A. Filter strawberry slurry through cheesecloth ____ To precipitate DNA from solution
- B. Mash strawberry with salty/soapy solution ____ Separate components of the cell
- C. Initial smashing and grinding of strawberry ____ Break open the cells
- D. Addition of ethanol to filtered extract ____ Break up proteins and dissolve cell Membranes

2. What did the DNA look like? Relate what you know about the chemical structure of DNA to what you observed today.

3. Explain what happened in the final step when you added ethanol to your strawberry extract. (*Hint: DNA is soluble in water, but not in ethanol*)

4. A person cannot see a single cotton thread 100 feet away, but if you wound thousands of threads together into a rope, it would be visible much further away. Is this statement analogous to our DNA extraction? Explain.

5. Why is it important for scientists to be able to remove DNA from an organism? List two reasons.

6. Is there DNA in your food? _____ How do you know?